

AMPRIA INSTRUCTIONS

This device has been provided to you to protect your amputation site as well as to keep your knee straight.



This device should be high enough on the leg so that ½ of the white cup at the end extends below the edge of the blue sleeve.

To take off the device, open only the front bottom flap and tabs as shown in the picture with red "X".



Your kneecap should be in the same area as the white dot.

Supracondylar strap sits ABOVE the kneecap for soft suspension. This is necessary to prevent the device from falling down.

Stiff plastic on the inside and outside of the leg to reduce knee bending.

Keep clover leaf tabs stuck to device on the inside, outside, and back (not shown). Only un-velcro the top one marked by the "X".

White foam pad shows BELOW the blue neoprene edge. This protects the end of the limb in case of a fall.

Please call Optimus Prosthetics with any questions.

Dayton (937) 454-1900 * Columbus (614) 263-5462 * Cincinnati (513) 918-2320